

Anti-hnRNP K Rabbit antibody



Catalog Number: 168939

Product name

Anti-hnRNP K Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

hnRNP K Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human hnRNP K

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions

WB: 1/1000-1/5000

IHC: 1/20-1/50

Validations



Western blot detection of hnRNP K in K562,PC3,HeLa cell lysates using hnRNP K Rabbit mAb(1:500 diluted).Predicted band size:51kDa.Observed band size:62kDa.