

EMILIN1 antibody



Catalog Number: 110239

Product name

EMILIN1 antibody

Specificity

human; other species not tested.

Antibody description

EMILIN1 Rabbit Polyclonal antibody. Positive IHC detected in human ovary tumor. Positive WB detected in human kidney tissue, human brain tissue, human placenta tissue. Observed molecular weight by Western-blot: 120kd

Preparation

This antibody was obtained by immunization of EMILIN1 recombinant protein (Accession Number: BC007530). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



human kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:110239(EMILIN1 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:110239(EMILIN1 antibody) at dilution of 1:100 (under 10x lens)