

HMGA2 Polyclonal Antibody



Catalog Number: 161239

Product name

HMGA2 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to HMGA2

Preparation

Antigen: A synthetic peptide of human HMGA2

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, FC

Dilutions

WB 1:500 - 1:2000

FC 1:20 - 1:50

Validations



Western blot - HMGA2 Polyclonal Antibody

Western blot analysis of extracts of HepG2 cells, using HMGA2 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.