

Anti-LGALS3 antibody



Catalog Number: 106239

Product name

Anti-LGALS3 antibody

Specificity

Human Galectin-3 / LGALS3

No cross-reactivity in ELISA with Human Galectin-1 / LGALS1; Human Galectin-7 / LGALS7; Human Galectin-8 / LGALS8

Antibody description

Rabbit monoclonal to LGALS3

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human Galectin-3 / LGALS3 (rh Galectin-3 / LGALS3; P17931; Met1-Ile250).

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions


IHC-P: 1-10 µg/mL

Validations



Galectin-3 / LGALS3 Antibody, Rabbit MAb, Immunocytochemistry

Immunocytochemical staining of human Galectin3 in human lung with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded



sections).

Galectin-3 / LGALS3 Antibody, Rabbit MAb, Immunocytochemistry

Immunocytochemical staining of human Galectin3 in human colon carcinoma with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections).