

Desmin antibody

Catalog Number: 107239

Product name

Desmin antibody

WB: 1:500-1:5000

IHC: 1:50-1:200

Specificity

Human, Mouse; other species not tested.

Antibody description

Desmin Mouse Monoclonal antibody. Positive IHC detected in human heart tissue, human hysteryomyoma tissue, human skeletal muscle tissue. Positive WB detected in Human heart. Observed molecular weight by Western-blot: 55 kDa

Preparation

This antibody was obtained by immunization of Desmin recombinant protein (Accession Number: NM_001927). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a


Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:107239(DES antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human heart using Catalog No:107239(DES antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human heart using Catalog No:107239(DES antibody) at dilution of 1:50 (under 40x lens)