

Mouse Tnf/Tnfa/Tnfsf2 recombinant protein



Catalog Number: 518239

General Information

Gene Name Synonym

Tumor Necrosis Factor; Cachectin; TNF-Alpha;
Tumor Necrosis Factor Ligand Superfamily
Member 2; TNF-a; Tumor Necrosis Factor,
Membrane Form; Tumor Necrosis Factor, Soluble
Form; Tnf; Tnfa; Tnfsf2

Protein Construction

A DNA sequence encoding the Mouse
Tnf/Tnfa/Tnfsf2 Asp89-Leu235 is expressed.

Organism

Mouse

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL
method.

Stability

Store for up to 12 months at -20°C to -80°C as
lyophilized powder. Aliquots of reconstituted
samples are stable at < -20°C for 3 months.

Formulation

Lyophilized from a 0.2 µm filtered solution of PBS,
pH 7.4.

Usage Guide


Storage

Lyophilized protein should be stored at <-20°C,
though stable at room temperature for 3 weeks.
Reconstituted protein solution can be stored at
4-7°C for 2-7 days.

Reconstitution

Always centrifuge tubes before opening. Do not
mix by vortex or pipetting. It is not recommended
to reconstitute to a concentration less than
100µg/ml. Dissolve the lyophilized protein in
ddH₂O. Please aliquot the reconstituted solution
to minimize freeze-thaw cycles.

SDS-PAGE



Mouse Tnf/Tnfa/Tnfsf2 recombinant protein