

# DHCR7 Polyclonal Antibody



Catalog Number: 160639

## Product name

DHCR7 Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to DHCR7

## Preparation

Antigen: A synthetic peptide of human DHCR7

## Formulation

PBS with 0.02% sodium azide, pH7.3.

## Storage

Store at 4°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100

IF 1:50 - 1:100

## Validations



## Western blot - DHCR7 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using DHCR7 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.