

# MEFV antibody



Catalog Number: 112639

## Product name

MEFV antibody

## Specificity

Human; other species not tested.

## Antibody description

MEFV Rabbit Polyclonal antibody. Positive WB detected in human plasma tissue. Observed molecular weight by Western-blot: 86 kDa

## Preparation

This antibody was obtained by immunization of MEFV recombinant protein (Accession Number: NM\_000243). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human plasma tissue were subjected to SDS PAGE followed by western blot with Catalog No:112639(MEFV Antibody) at dilution of 1:1000