

# TAPBP Polyclonal Antibody



Catalog Number: 163639

## Product name

TAPBP Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to TAPBP

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 30-260 of human TAPBP (NP\_757345.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100


IF 1:50 - 1:100

## Validations




Western blot - TAPBP Polyclonal Antibody

Western blot analysis of extracts of HT-29 cells, using TAPBP antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - TAPBP Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using TAPBP Antibody at dilution



of 1:200 (40x lens).

# TAPBP Polyclonal Antibody




Catalog Number: 163639

---

Immunohistochemistry - TAPBP Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat testis using TAPBP Antibody at dilution of 1:200



(40x lens).

Immunohistochemistry - TAPBP Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human gastric cancer using TAPBP Antibody at dilution of 1:200 (40x lens).



Immunohistochemistry - TAPBP Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse testis using TAPBP Antibody at dilution of 1:200 (40x lens).