

Anti-IGLC1 antibody



Catalog Number: 176639

Product name

Anti-IGLC1 antibody

Specificity

Human

Antibody description

Mouse monoclonal antibody to IGLC1

Preparation

This antigen of this antibody was native protein.

Formulation

Liquid, 1*TBS (pH7.4), 0.2% BSA, 50% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications

WB, ELISA

Dilutions

WB: 1:1000

ELISA: Use at an assay dependent concentration.

Validations




Fig1.; Western blot analysis of Human Iambda light chain on human IgG protein lysate. Proteins were transferred to a PVDF membrane and blocked with 5% BSA in PBS for 1 hour at room temperature. The primary antibody (1/1000) was used in 5% BSA at room temperature for 2 hours. Goat Anti-Mouse IgG - HRP Secondary Antibody (HA1006) at 1:5,000 dilution was used for 1 hour at room temperature.