

Human UCMA/C10orf49 recombinant protein



Catalog Number: 519639

General Information

Gene Name Synonym

Unique cartilage matrix-associated protein, UCMA, Gla-rich protein, C10orf49

Protein Construction

A DNA sequence encoding the Human UCMA/C10orf49 Ser65-Thr138 is expressed.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM PB \square 150mM NaCl, pH7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human UCMA/C10orf49 recombinant protein