

Eotaxin antibody

Catalog Number: 110339

Product name

Eotaxin antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Eotaxin Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue, mouse lung tissue, mouse small intestine tissue, mouse thymus tissue. Observed molecular weight by Western-blot: 10kd, 20kd

Preparation

This antibody was obtained by immunization of Eotaxin recombinant protein (Accession Number: NM_002986). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:110339(Eotaxin antibody) at dilution of 1:1000