

LGR6 Polyclonal Antibody



Catalog Number: 165339

Product name

LGR6 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to LGR6

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 350-510 of human LGR6 (NP_067649.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations



Western blot - LGR6 Polyclonal Antibody


Western blot analysis of extracts of various cell lines, using LGR6 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure



time: 90s.

Immunohistochemistry - LGR6 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using LGR6 Antibody at dilution of 1:100



(40x lens).

LGR6 Polyclonal Antibody



Catalog Number: 165339

Immunohistochemistry - LGR6 Polyclonal Antibody mouse brain using LGR6 Antibody at dilution of 1:100 (40x lens).
Immunohistochemistry of paraffin-embedded