

# Anti-ASGPR1 antibody



Catalog Number: 106339

## Product name

Anti-ASGPR1 antibody

## Immunogen

[Mouse ASGPR1 \(His Tag\) recombinant protein](#)

## Specificity

Mouse ASGPR1 / ASGR1

## Antibody description

Rabbit monoclonal to ASGPR1

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse ASGPR1 / ASGR1 (rM ASGPR1 / ASGR1; NP\_033844.1; Ser60-Asn284).

## Formulation

0.2 µm filtered solution in PBS with 5% trehalose

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions

IHC-P: 0.5-2 µg/mL

## Validations



ASGPR1 / ASGR1 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of mouse ASGR1 in mouse liver with rabbit monoclonal antibody (2 µg/mL, formalin-fixed paraffin embedded sections).