

Anti-TCFL5 antibody



Catalog Number: 176339

Product name

Anti-TCFL5 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to TCFL5

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tcfl5. aa range:311-360

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1:: Western Blot analysis of RAW264.7 cells using TCFL5 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventibiotech, MN, USA).