

# PCK1 antibody

Catalog Number: 113739

## Product name

PCK1 antibody

IP: 1:200-1:2000

IF: 1:20-1:200

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PCK1 Rabbit Polyclonal antibody. Positive IP detected in HepG2 cells. Positive WB detected in HepG2 cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 67 kDa

## Preparation

This antibody was obtained by immunization of PCK1 recombinant protein (Accession Number: NM\_002591). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IP, IF

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:113739(PCK1 antibody) at dilution of 1:1000



IP Result of anti-PCK1 (IP:Catalog No:113739, 4ug; Detection:Catalog No:113739 1:500) with HepG2 cells lysate 3600ug.



Immunofluorescent analysis of (10% Formaldehyde) fixed HepG2 cells using Catalog

# PCK1 antibody



Catalog Number: 113739

---

No:113739(PCK1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)