

# TRPM2 Polyclonal Antibody



Catalog Number: 164739

## Product name

TRPM2 Polyclonal Antibody

## Specificity

Mouse, Rat

## Antibody description

Polyclonal antibody to TRPM2

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 900-1000 of human TRPM2 (NP\_003298.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IF

## Dilutions


IF 1:50 - 1:100

## Validations



Immunofluorescence - TRPM2 Polyclonal Antibody

Immunofluorescence analysis of rat brain using TRPM2 antibody at dilution of 1:100. Blue: DAPI



for nuclear staining.

Immunofluorescence - TRPM2 Polyclonal Antibody

Immunofluorescence analysis of mouse brain using TRPM2 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.