

Anti-DNAJB13 antibody



Catalog Number: 175739

Product name

Anti-DNAJB13 antibody

Specificity

Human, Mouse, Rat, Dog, Pig, Cow, Horse, Rabbit, Guinea Pig

Antibody description

Rabbit polyclonal antibody to DNAJB13

Preparation

This antigen of this antibody was klh conjugated synthetic peptide derived from human tsarg6 185-280/316

Formulation

Liquid, 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

Storage

Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB:1:500-2000

Validations





Fig1: Sample:; U2OS(Human) Cell Lysate at 30 ug; HeLa(Human) Cell Lysate at 30 ug; Primary: Anti- TSARG6 at 1/1000 dilution; Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution; Predicted band size: 36 kD; Observed



band size: 34 kD

Fig2: Sample:; U87MG(Human) Cell Lysate at 30 ug; Primary: Anti- TSARG6 at 1/1000 dilution; Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution; Predicted band size: 36 kD; Observed band size: 35 kD