

Anti-ACTR1B Rabbit antibody



Catalog Number: 167739

Product name

Anti-ACTR1B Rabbit antibody

Specificity

Human

Antibody description

ACTR1B Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human ACTR1B

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/2000


IHC: 1/50

ICC/IF: 1/50

FC: 1/50

IP: 1/20

Validations



Western blot detection of ACTR1B in HeLa, A549, HL-60, U251, U87-MG cell lysates using ACTR1B Rabbit mAb(1:1000 diluted). Predicted band size: 42kDa. Observed band size: 42kDa.