

NCX1 antibody

Catalog Number: 113049

Product name

NCX1 antibody

Specificity

Human, Rat, Mouse; other species not tested.

Antibody description

NCX1 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue, rat brain tissue. Observed molecular weight by Western-blot: 120 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_021097). Purification method: Antigen affinity purified.

Formulation

0.1M NaHCO₃, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113049(SLC8A1 Antibody) at dilution of



1:1000

IP Result of anti-SLC8A1 (IP:Catalog No:113049, 4ug; Detection:Catalog No:113049 1:1000) with mouse brain tissue lysate 4000ug.