

Anti-LAMA5 antibody

Catalog Number: 176049

Product name

Anti-LAMA5 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to LAMA5

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human lama5. aa range:2381-2430

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.: Western Blot analysis of RAW HepG2 A549 cells using Laminin α -5 Polyclonal Antibody