

# Marcks antibody

Catalog Number: 112449

## Product name

Marcks antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

Marcks Rabbit Polyclonal antibody. Positive IHC detected in human prostate cancer tissue, mouse pancreas tissue. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 88kd

## Preparation

This antibody was obtained by immunization of Recombinant Protein (Accession Number: NM\_008538). Purification method: Protein A purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112449(Marcks antibody) at dilution of 1:2000



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:112449(Marcks antibody) at dilution of 1:200 (under 10x lens)