

# NVL antibody



Catalog Number: 113449

## Product name

NVL antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

NVL Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in mouse lung tissue, mouse kidney tissue, mouse thymus tissue. Observed molecular weight by Western-blot: 95kd

## Preparation

This antibody was obtained by immunization of NVL recombinant protein (Accession Number: BC012105). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:113449(NVL antibody) at dilution of 1:600



Immunofluorescent analysis of HepG2 cells using Catalog No:113449(NVL Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG