

Phospho-REL-S503 pAb



Catalog Number: 166449

Product name

Phospho-REL-S503 pAb

Specificity

Human

Antibody description

Polyclonal antibody to Phospho-REL-S503

Preparation

Antigen: A phospho specific peptide corresponding to residues surrounding S503 of human REL

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:100 - 1:200

Validations



Immunofluorescence - Phospho-REL-S503 pAb

Immunofluorescence analysis of methanol-fixed HeLa cells using Phospho-REL-S503 antibody .