

# Anti-Ku70 (6H10) Mouse antibody



Catalog Number: 168449

## Product name

Anti-Ku70 (6H10) Mouse antibody

## Dilutions

WB: 1/1000

## Specificity

;ICC: 1/200

Human

## Validations

## Antibody description

Ku70 (6H10) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human Ku70 protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality


Monoclonal

## Ig Type


Mouse IgG2b

## Applications

WB, ICC/IF



Immunocytochemistry staining of HeLa cells fixed with 4% Paraformaldehyde and using Ku70 mouse



mAb (dilution 1:200).

Western blot detection of Ku70 in HeLa, A549, MCF7, Jurkat and Raji cell lysates using Ku70 mouse mAb (1:1000 diluted). Predicted band size: 70KDa. Observed band size: 70KDa.