

GALNT9 antibody



Catalog Number: 110849

Product name

GALNT9 antibody

Immunogen

[Human GALNT9 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GALNT9 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, C6 cells. Observed molecular weight by Western-blot: 66-68 kDa

Preparation

This antibody was obtained by immunization of GALNT9 recombinant protein (Accession Number: BC157845). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110849(GALNT9 Antibody) at dilution of 1:1000