

# GALNT9 antibody



Catalog Number: 110849

## Product name

GALNT9 antibody

## Immunogen

[Human GALNT9 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

GALNT9 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, C6 cells. Observed molecular weight by Western-blot: 66-68 kDa

## Preparation

This antibody was obtained by immunization of GALNT9 recombinant protein (Accession Number: BC157845). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110849(GALNT9 Antibody) at dilution of 1:1000