

Anti-C2CD6 antibody



Catalog Number: 175849

Product name

Anti-C2CD6 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to C2CD6

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human als2cr11. aa range:491-540

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of 3T3 cells using ALS2CR11 Polyclonal Antibody