

Anti-Angiopoietin 1 Rabbit antibody



Catalog Number: 168849

Product name

Anti-Angiopoietin 1 Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Rat

Applications

WB

Antibody description

Angiopoietin 1 Rabbit monoclonal antibody

Dilutions

1:2000-1:10000


Preparation

Antigen: A synthetic peptide of human Angiopoietin 1

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of Angiopoietin 1 in Jurkat, C6, HeLa cell lysates using Angiopoietin 1 Rabbit mAb(1:1000 diluted). Predicted band size: 58kDa. Observed band size: 75kDa.

Clonality

Monoclonal