

Rat Apolipoprotein H/APOH (His Tag)



Catalog Number: 500149

General Information

Gene Name Synonym

Apolipoprotein H (Beta-2-glycoprotein I); Beta-2-glycoprotein 1

Protein Construction

A DNA sequence encoding the rat Apoh (NP_001009626.1) (Met1-Cys345) was expressed with a polyhistidine tag at the C-terminus.

Organism

Rat

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE.

Endotoxin

< 1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gly 20

Molecular Mass

The recombinant rat Apoh consists 337 amino acids and predicts a molecular mass of 37.9 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Rat Apolipoprotein H / APOH Protein (His Tag)