

# SFRS2 antibody

Catalog Number: 115149

## Product name

SFRS2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SFRS2 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue, A549 cells, mouse testis tissue, rat brain tissue. Observed molecular weight by Western-blot: 25-35 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_003016). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP

## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115149(SFRS2 antibody) at dilution of 1:300



IP Result of anti-SFRS2 (IP:Catalog No:115149, 4ug; Detection:Catalog No:115149 1:300) with mouse brain tissue lysate 4000ug.