

Anti-OR10H4 antibody



Catalog Number: 176149

Product name

Anti-OR10H4 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to OR10H4

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human or10h4. aa range:161-210

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1.; Western Blot analysis of RAW264.7 cells using Olfactory receptor 10H4 Polyclonal Antibody