

TNNC2 antibody

Catalog Number: 116149

Product name

TNNC2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TNNC2 Rabbit Polyclonal antibody. Positive IHC detected in human skeletal muscle tissue. Positive WB detected in mouse skeletal muscle tissue, rat skeletal muscle tissue. Observed molecular weight by Western-blot: 20-23 kDa

Preparation

This antibody was obtained by immunization of TNNC2 recombinant protein (Accession Number: NM_003279). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:

WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:116149(TNNC2 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using Catalog No:116149(TNNC2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using Catalog

TNNC2 antibody



Catalog Number: 116149

No:116149(TNNC2 Antibody) at dilution of 1:50 (under 40x lens)