

HSD17B7 antibody



Catalog Number: 111549

Product name

HSD17B7 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

HSD17B7 Rabbit Polyclonal antibody. Positive IP detected in MCF-7 cells. Positive WB detected in MCF7 cells, HEK-293 cells. Observed molecular weight by Western-blot: 35-38 kDa

Preparation

This antibody was obtained by immunization of HSD17B7 recombinant protein (Accession Number: NM_016371). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:111549(HSD17B7 antibody) at dilution of



IP Result of anti-HSD17B7 (IP:Catalog No:111549, 4ug; Detection:Catalog No:111549 1:1000) with MCF-7 cells lysate 3200ug.