

RAP2A antibody



Catalog Number: 114549

Product name

RAP2A antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

RAP2A Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, mouse brain tissue, PC-3 cells. Observed molecular weight by Western-blot: 21kd

Preparation

This antibody was obtained by immunization of RAP2A recombinant protein (Accession Number: NM_021033). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:114549(RAP2A antibody) at dilution of 1:1000