

NAPRT1 antibody

Catalog Number: 112949

Product name

NAPRT1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NAPRT1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HepG2 cells, human colon tissue, human small intestine tissue, L02 cells, mouse kidney tissue, mouse liver tissue. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 51 kDa, 56 kDa

Preparation

This antibody was obtained by immunization of NAPRT1 recombinant protein (Accession Number: NM_001363146). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:112949(NAPRT1 antibody) at dilution of



1:1000

IP Result of anti-NAPRT1 (IP:Catalog No:112949, 4ug; Detection:Catalog No:112949 1:1000) with HEK-293 cells lysate 2800ug.