

# Anti-Ferritin Heavy Chain Rabbit antibody



Catalog Number: 166949

## Product name

Anti-Ferritin Heavy Chain Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Ferritin Heavy Chain Rabbit polyclonal antibody

## Preparation

Antigen: A synthesized peptide derived from human Ferritin

## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, ICC/IF

## Dilutions


WB: 1/500-1/2000

ICC/IF: 1/50-1/200

## Validations



Western blot analysis of Ferritin expression in



Jurkat cell lysate.

Immunofluorescent analysis of Jurkat cells, using Ferritin Antibody.