

ZNF143 antibody



Catalog Number: 116949

Product name

ZNF143 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ZNF143 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, mouse brain tissue. Observed molecular weight by Western-blot: 68 kDa

Preparation

This antibody was obtained by immunization of ZNF143 recombinant protein (Accession Number: BC020219). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:116949(ZNF143 antibody) at dilution of 1:800