

Anti-SBDS Rabbit antibody



Catalog Number: 167949

Product name

Anti-SBDS Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

SBDS Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human SBDS

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, IP


Dilutions

WB: 1/2000-1/10000

ICC/IF: 1/50

IP: 1/20

Validations



Western blot detection of SBDS in A549,HL-60,U2OS,U87-MG cell lysates using SBDS Rabbit mAb(1:1000 diluted).Predicted band size:29kDa.Observed band size:29kDa.