

Anti-Phospholipase C beta 1/PLCB1 Rabbit antibody



Catalog Number: 168949

Product name

Anti-Phospholipase C beta 1/PLCB1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Phospholipase C beta 1/PLCB1 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of mouse Phospholipase C beta 1

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions

WB: 1/1000

IHC: 1/200

Validations



Western blot detection of Phospholipase C beta 1 in Rat Brain lysates using Phospholipase C beta 1 Rabbit mAb(1:500 diluted). Predicted band size:138kDa. Observed band size:138kDa.