

PNRC2 Polyclonal Antibody



Catalog Number: 163249

Product name

PNRC2 Polyclonal Antibody

Specificity

Mouse

Antibody description

Polyclonal antibody to PNRC2

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PNRC2 (NP_060231.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - PNRC2 Polyclonal Antibody

Western blot analysis of extracts of mouse lung, using PNRC2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 3min.