

Anti-GRB2 (5G3) Mouse antibody



Catalog Number: 168249

Product name

Anti-GRB2 (5G3) Mouse antibody

Ig Type

Mouse IgG1

Specificity

Human, Mouse, Rat

Applications

WB

Antibody description

GRB2 (5G3) Mouse monoclonal antibody

Dilutions

WB: 1/1000

Preparation

Antigen: Recombinant protein of human GRB2.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.


Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Validations



Western blot detection of GRB2 in K562, Raji, Jurkat, Ramos, C6, 3T3 and HeLa cell lysates using GRB2 mouse mAb (1:1000 diluted). Predicted band size: 25KDa. Observed band size: 25KDa.