

# Anti-GRB2 (5G3) Mouse antibody



Catalog Number: 168249

## Product name

Anti-GRB2 (5G3) Mouse antibody

## Specificity

Human, Mouse, Rat

## Antibody description

GRB2 (5G3) Mouse monoclonal antibody

## Preparation

Antigen: Recombinant protein of human GRB2.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of GRB2 in K562, Raji, Jurkat, Ramos, C6, 3T3 and HeLa cell lysates using GRB2 mouse mAb (1:1000 diluted). Predicted band size: 25KDa. Observed band size: 25KDa.