

Anti-Daxx Rabbit antibody



Catalog Number: 167649

Product name

Anti-Daxx Rabbit antibody

Specificity

Human

Antibody description

Daxx Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Daxx

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


WB: 1/2000

IHC: 1/20

ICC/IF: 1/20

FC: 1/20

Validations



Western blot detection of Daxx in
HepG2,U251,HL-60 using Daxx Rabbit
mAb(1:1000 diluted)