

# LSM4 antibody



Catalog Number: 112349

## Product name

LSM4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

LSM4 Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue. Positive WB detected in MCF7 cells, human placenta tissue, MDA-MB-231 cells. Observed molecular weight by Western-blot: 15 kDa

## Preparation

This antibody was obtained by immunization of LSM4 recombinant protein (Accession Number: NM\_012321). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:112349(LSM4 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:112349(LSM4 antibody) at dilution of 1:100 (under 10x lens)