

SLC7A5 antibody



Catalog Number: 115349

Product name

SLC7A5 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SLC7A5 Rabbit Polyclonal antibody. Positive IHC detected in human lung tissue, human stomach cancer tissue, human testis tissue. Positive WB detected in mouse lung tissue. Observed molecular weight by Western-blot: 45-50 kDa

Preparation

This antibody was obtained by immunization of SLC7A5 recombinant protein (Accession Number: NM_003486). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human lung using Catalog No:115349(SLC7A5 antibody) at dilution of 1:50 (under 10x lens)



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:115349(SLC7A5 Antibody) at dilution of 1:600