

CLEC16A antibody



Catalog Number: 109349

Product name

CLEC16A antibody

Immunogen

[Human CLEC16A Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

CLEC16A Rabbit Polyclonal antibody. Positive WB detected in mouse spleen tissue. Observed molecular weight by Western-blot: 104-110 kDa

Preparation

This antibody was obtained by immunization of CLEC16A recombinant protein (Accession Number: NM_001243403). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:109349(CLEC16A Antibody) at dilution of 1:600