

IL6R antibody

Catalog Number: 111749

Product name

IL6R antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

IL6R Rabbit Polyclonal antibody. Positive WB detected in mouse spleen tissue. Positive IP detected in mouse spleen tissue. Observed molecular weight by Western-blot: 80 kDa

Preparation

This antibody was obtained by immunization of IL6R recombinant protein (Accession Number: NM_000565). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:111749(IL6R antibody) at dilution of 1:1000



IP Result of anti-IL6R (IP:Catalog No:111749, 4ug; Detection:Catalog No:111749 1:500) with mouse spleen tissue lysate 4000ug.