

# MPP5 antibody

Catalog Number: 112749

## Product name

MPP5 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

MPP5 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in Y79 cells, human brain tissue. Observed molecular weight by Western-blot: 70-75kd

## Preparation

This antibody was obtained by immunization of MPP5 recombinant protein (Accession Number: XM\_047431691). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:112749(MPP5 antibody) at dilution of 1:400



IP Result of anti-PALS1 (IP:Catalog No:112749, 3ug; Detection:Catalog No:112749 1:1000) with mouse brain tissue lysate 8000ug.