

# VHL antibody



Catalog Number: 116749

---

## Product name

VHL antibody

## Specificity

Human; other species not tested.

## Antibody description

VHL Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells.

## Preparation

This antibody was obtained by immunization of VHL recombinant protein (Accession Number: NM\_198156). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IF

## Dilutions

Recommended Dilution:

IF: 1:10-1:100

## Validations



Immunofluorescent analysis of HepG2 cells using Catalog No:116749(VHL Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)